



INSTEP LASER CANE



SKU: US-IN-LC

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

lightweight height adjustable aluminium cane that projects a battery operated bright red line across the user's path. The laser beam is an obstacle free visual cue, which can enable people with Parkinson's disease or similar conditions, to break "freezing" episodes and increase stride length. The laser light appears when the user presses down on the cane and disappears when they lift it. The laser is powered by two AA batteries.

FEATURES

- For Parkinson's freezing and other neurological conditions
- Projects a bright red line across your path as a safe, obstacle-free visual cue
- Enables Parkinson's patients and others to break "freezing" episodes and increase stride length
- Lightweight
- Height adjustable
- Powered by AA batteries

SPECIFICATIONS

- Handle height 790 - 1016mm
- Height adjustments in 25mm increments

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au