



INVACARE COLIBRI MOBILITY SCOOTER

~~\$2,297.00~~ Original price was:

~~\$2,297.00.~~ ~~\$2,067.04~~ Current price is: \$2,067.04. Sale

Available Online Only!



SKU: 1577246-JR-H

PRODUCT DESCRIPTION

Purpose & Benefits of Invacare Colibri Mobility Scooter

The **Invacare Colibri Mobility Scooter** is a stylish and compact 4-wheel scooter designed for users who need a lightweight, easy-to-transport mobility solution. Ideal for everyday use in indoor settings or smooth outdoor surfaces, the Colibri offers excellent manoeuvrability, user comfort, and quick disassembly with Invacare's patented LiteLock™ system.

This scooter is perfect for travel, shopping, or social outings, offering a comfortable padded seat, intuitive controls, and tool-free assembly and disassembly. Its compact size and portability make it especially suitable for users seeking independence and convenience on the go.

Key Features

- Compact and lightweight design ideal for transport and storage
- Tool-free LiteLock™ system for fast disassembly
- Comfortable padded seat with swivel function
- Ergonomic tiller with intuitive controls
- Suitable for indoor use and smooth outdoor surfaces

- Easily fits into most car boots
- Designed for portability, independence, and daily convenience

Specifications

- Weight capacity: 136 kg
- Product weight: 49 kg
- Range up to 16 km
- Max speed: 8 km/h
- Turning Radius: 2200 mm
- Batteries: 2 x 18 Ah
- Width: 465 mm
- Length: 1010 mm