



## INVACARE COLIBRI MOBILITY SCOOTER



**SKU:** 1577246-JR-H

## PRODUCT DESCRIPTION

### Purpose & Benefits of Invacare Colibri Mobility Scooter

The **Invacare Colibri Mobility Scooter** is a stylish and compact 4-wheel scooter designed for users who need a lightweight, easy-to-transport mobility solution. Ideal for everyday use in indoor settings or smooth outdoor surfaces, the Colibri offers excellent manoeuvrability, user comfort, and quick disassembly with Invacare's patented LiteLock™ system.

This scooter is perfect for travel, shopping, or social outings, offering a comfortable padded seat, intuitive controls, and tool-free assembly and disassembly. Its compact size and portability make it especially suitable for users seeking independence and convenience on the go.

### Key Features

- Compact and lightweight design ideal for transport and storage
- Tool-free LiteLock™ system for fast disassembly
- Comfortable padded seat with swivel function
- Ergonomic tiller with intuitive controls
- Suitable for indoor use and smooth outdoor surfaces

- Easily fits into most car boots
- Designed for portability, independence, and daily convenience

## Specifications

- Weight capacity: 136 kg
- Product weight: 49 kg
- Range up to 16 km
- Max speed: 8 km/h
- Turning Radius: 2200 mm
- Batteries: 2 x 18 Ah
- Width: 465 mm
- Length: 1010 mm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)

