



INVACARE COMET ALPINE PLUS SCOOTER

\$6,898.25



SKU: 1637602-B-H

PRODUCT DESCRIPTION

Purpose & Benefits of Invacare Comet Alpine Plus Scooter

The **Invacare Comet Alpine Plus Scooter** is a high-performance, all-terrain mobility scooter designed for users who require stability, power, and confidence over rough or uneven ground. Built for outdoor exploration, the Comet Alpine Plus features large wheels, a high ground clearance, and robust suspension, making it ideal for rural areas, hilly paths, and challenging environments.

With a high-capacity motor and battery, this scooter provides extended range and reliable performance, while its ergonomic seating and intelligent safety systems ensure user comfort and peace of mind. Whether navigating slopes or cruising through the countryside, the Comet Alpine Plus offers freedom and independence without compromise.

Key Features

- Heavy-duty design for rough terrain and outdoor conditions
- Large pneumatic tyres and high ground clearance for superior traction
- Powerful motor and battery for extended range and hill climbing
- Full suspension system for a smooth and stable ride
- High-back adjustable seat for long-term comfort

- Ergonomic controls with safety features including automatic cornering speed reduction
- Ideal for rural, hilly, or uneven environments

Specifications

- **Weight Capacity:** 150kg
- **Weight:** 143kg
- **Width:** 66.5-68.5cm
- **Length:** 145cm
- **Backrest Height:** 63cm
- **Backrest Width:** 51cm
- **Seat Depth:** 46cm
- **Floor-to-Seat Height:** 44-51cm
- **Ground Clearance:** 10cm
- **Max Climb Angle:** 10 degrees
- **Speed:** 10km/h
- **Range:** Up to 58km