




INVACARE COMET ALPINE PLUS SCOOTER



SKU: 1637602-B-H

PRODUCT DESCRIPTION

Purpose & Benefits of Invacare Comet Alpine Plus Scooter

The **Invacare Comet Alpine Plus Scooter** is a high-performance, all-terrain mobility scooter designed for users who require stability, power, and confidence over rough or uneven ground. Built for outdoor exploration, the Comet Alpine Plus features large wheels, a high ground clearance, and robust suspension, making it ideal for rural areas, hilly paths, and challenging environments.

With a high-capacity motor and battery, this scooter provides extended range and reliable performance, while its ergonomic seating and intelligent safety systems ensure user comfort and peace of mind. Whether navigating slopes or cruising through the countryside, the Comet Alpine Plus offers freedom and independence without compromise.

Key Features

- Heavy-duty design for rough terrain and outdoor conditions
- Large pneumatic tyres and high ground clearance for superior traction
- Powerful motor and battery for extended range and hill climbing
- Full suspension system for a smooth and stable ride
- High-back adjustable seat for long-term comfort

- Ergonomic controls with safety features including automatic cornering speed reduction
- Ideal for rural, hilly, or uneven environments

Specifications

- **Weight Capacity:** 150kg
- **Weight:** 143kg
- **Width:** 66.5-68.5cm
- **Length:** 145cm
- **Backrest Height:** 63cm
- **Backrest Width:** 51cm
- **Seat Depth:** 46cm
- **Floor-to-Seat Height:** 44-51cm
- **Ground Clearance:** 10cm
- **Max Climb Angle:** 10 degrees
- **Speed:** 10km/h
- **Range:** Up to 58km

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au