




INVACARE FOLDING SIDE RAILS (PAIR) - CS8 BED COMPATIBLE



SKU: N/A

PRODUCT DESCRIPTION

Purpose & Benefits

The **Invacare Folding Side Rails (Pair)** are specifically designed to provide enhanced safety and fall prevention for users of the **Invacare CS8 Floorline Care Bed**. Ideal for both home and clinical environments, these bolt-on $\frac{3}{4}$ length rails are easy to operate and fold down for carer access or patient transfers.

Built for durability and ease of use, the folding rails integrate seamlessly with the CS8 bed frame, supporting patient reassurance and mobility without compromising functionality or comfort.

Features

- $\frac{3}{4}$ length folding side rails for the Invacare CS8 bed
- Helps prevent falls and supports repositioning
- Smooth folding action allows quick access and transfers
- Secure bolt-on design for stable attachment
- Designed for use in aged care, hospital, and home care environments
- Easy to clean with minimal maintenance

Specifications

- Compatibility: Invacare CS8 Electrically Actuated Care Bed
- Mounting: Bolt-on attachment
- Configuration: Pair (left and right rails)
- Fold Type: Swing-down mechanism
- Length: $\frac{3}{4}$ bed length coverage
- Material: Powder-coated steel and durable plastic components
- Cleaning: Wipe-clean surfaces
- Manufacturer: Invacare

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au