




INVACARE FOLDING SIDE RAILS (PAIR) - CS8 BED COMPATIBLE



SKU: N/A

PRODUCT DESCRIPTION

Purpose & Benefits

The **Invacare Folding Side Rails (Pair)** are specifically designed to provide enhanced safety and fall prevention for users of the **Invacare CS8 Floorline Care Bed**. Ideal for both home and clinical environments, these bolt-on $\frac{3}{4}$ length rails are easy to operate and fold down for carer access or patient transfers.

Built for durability and ease of use, the folding rails integrate seamlessly with the CS8 bed frame, supporting patient reassurance and mobility without compromising functionality or comfort.

Features

- $\frac{3}{4}$ length folding side rails for the Invacare CS8 bed
- Helps prevent falls and supports repositioning
- Smooth folding action allows quick access and transfers
- Secure bolt-on design for stable attachment
- Designed for use in aged care, hospital, and home care environments
- Easy to clean with minimal maintenance

Specifications

- Compatibility: Invacare CS8 Electrically Actuated Care Bed
- Mounting: Bolt-on attachment
- Configuration: Pair (left and right rails)
- Fold Type: Swing-down mechanism
- Length: $\frac{3}{4}$ bed length coverage
- Material: Powder-coated steel and durable plastic components
- Cleaning: Wipe-clean surfaces
- Manufacturer: Invacare