



INVACARE OCTAVE BARIATRIC BED WITH OVERHEAD POLE & SIDE RAILS



SKU: OCTAVE_C2AB1_M1_1-WM

PRODUCT DESCRIPTION

Purpose & Benefits

The **Invacare Octave Bariatric Bed** with overhead pole and side rails is a fully electric, four-section profiling bed designed to deliver maximum comfort, safety and independence for bariatric users. Engineered with strength, accessibility and ergonomic operation in mind, it provides smooth height and position adjustments while ensuring carers can assist patients safely and efficiently. Ideal for hospitals, aged care and home environments, the Octave enhances patient mobility and overall care quality.

Features

- Fully electric 4-section mattress platform for adjustable comfort and positioning
- Extended bed frame design offers full support across the width and edges of the bed base
- Detachable bed ends for quick access to both head and foot sections without removing side rails
- Extendable bed end allows flexible length adjustment (10-20 cm extension)
- Adjustable overhead lifting pole with handle that can be set to different heights and depths
- Integrated support handle assists users with safe and stable transfers in and out of bed
- Simple-to-use hand control for adjusting bed height and profiling functions

- Built-in transport rack for convenient relocation and compact storage
- Includes side rails for added patient security and safety

Specifications

Specification	Details
Model	Invacare Octave Bariatric Bed
Bed Type	Fully electric 4-section profiling bed
Included Accessories	Overhead lifting pole and side rails
Bed End Extension	Adjustable 10-20 cm extension range
Lifting Pole	Height and depth adjustable handle for user support
Frame Design	Reinforced steel frame with full-width support base
Control	Easy-to-use hand pendant for bed height and profiling functions
Transport System	Built-in transport rack for mobility and storage

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au