



INVACARE OCTAVE BARIATRIC BED WITH OVERHEAD POLE & SIDE RAILS



SKU: OCTAVE_C2AB1_M1_1-WM

PRODUCT DESCRIPTION

Purpose & Benefits

The **Invacare Octave Bariatric Bed** with overhead pole and side rails is a fully electric, four-section profiling bed designed to deliver maximum comfort, safety and independence for bariatric users. Engineered with strength, accessibility and ergonomic operation in mind, it provides smooth height and position adjustments while ensuring carers can assist patients safely and efficiently. Ideal for hospitals, aged care and home environments, the Octave enhances patient mobility and overall care quality.

Features

- Fully electric 4-section mattress platform for adjustable comfort and positioning
- Extended bed frame design offers full support across the width and edges of the bed base
- Detachable bed ends for quick access to both head and foot sections without removing side rails
- Extendable bed end allows flexible length adjustment (10-20 cm extension)
- Adjustable overhead lifting pole with handle that can be set to different heights and depths
- Integrated support handle assists users with safe and stable transfers in and out of bed
- Simple-to-use hand control for adjusting bed height and profiling functions



- Built-in transport rack for convenient relocation and compact storage
- Includes side rails for added patient security and safety

Specifications

Specification		Details
Model	Invacare Octave Bariatric Bed	

Bed Type Fully electric 4-section profiling bed Included Accessories Overhead lifting pole and side rails
Bed End Extension Adjustable 10–20 cm extension range

Lifting Pole Height and depth adjustable handle for user support Frame Design Reinforced steel frame with full-width support base

Control Easy-to-use hand pendant for bed height and profiling functions

Transport System Built-in transport rack for mobility and storage