



## INVACARE PEGASUS METRO SCOOTER



**SKU:** 1637600-JR-H

## PRODUCT DESCRIPTION

### Purpose & Benefits of Invacare Pegasus Metro Scooter

The **Invacare Pegasus Metro Scooter** is a versatile, mid-size mobility scooter designed for comfort, safety, and reliable performance during daily travel. Combining ergonomic features with intelligent electronics and a robust motor, the Pegasus Metro is perfect for users who need a balance of power and practicality for both urban and suburban environments.

This scooter offers excellent ride comfort with its adjustable suspension, padded high-back seat, and large pneumatic tyres. It also features Invacare's intelligent electronics system for smooth acceleration, automatic speed reduction on corners, and enhanced safety. Whether commuting locally or enjoying a longer outing, the Pegasus Metro delivers a confident and dependable driving experience.

### Key Features

- Mid-size design ideal for daily travel and moderate outdoor use
- Adjustable suspension and padded seat for enhanced ride comfort
- Large pneumatic tyres for smooth handling over uneven surfaces
- Automatic speed reduction for added cornering safety
- Ergonomic handlebars and intuitive controls

- Integrated front and rear lighting for increased visibility
- Reliable performance suited to both urban and suburban environments

## Specifications

<b>Total Length</b>	1270mm
<b>Driving Unit Width</b>	660mm
<b>Backrest Width</b>	510mm
<b>Seat Depth</b>	470mm
<b>Floor to Seat Height</b>	440mm - 510mm
<b>Backrest Height</b>	490mm
<b>Maximum Speed</b>	10kph
<b>Range</b>	54km
<b>Weight</b>	110kg
<b>SWL</b>	136kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)