



INVACARE PEGASUS METRO SCOOTER

\$5,188.21 Original price was: \$5,188.21.\$4,669.00Current price is: \$4,669.00.



SKU: 1637600-JR-H

PRODUCT DESCRIPTION

Purpose & Benefits of Invacare Pegasus Metro Scooter

The **Invacare Pegasus Metro Scooter** is a versatile, mid-size mobility scooter designed for comfort, safety, and reliable performance during daily travel. Combining ergonomic features with intelligent electronics and a robust motor, the Pegasus Metro is perfect for users who need a balance of power and practicality for both urban and suburban environments.

This scooter offers excellent ride comfort with its adjustable suspension, padded high-back seat, and large pneumatic tyres. It also features Invacare's intelligent electronics system for smooth acceleration, automatic speed reduction on corners, and enhanced safety. Whether commuting locally or enjoying a longer outing, the Pegasus Metro delivers a confident and dependable driving experience.

Key Features

- Mid-size design ideal for daily travel and moderate outdoor use
- Adjustable suspension and padded seat for enhanced ride comfort
- Large pneumatic tyres for smooth handling over uneven surfaces
- Automatic speed reduction for added cornering safety
- Ergonomic handlebars and intuitive controls



- Integrated front and rear lighting for increased visibility
- Reliable performance suited to both urban and suburban environments

Specifications

Total Length	1270mm	
Driving Unit Width	660mm	
Backrest Width	510mm	
Seat Depth	470mm	
Floor to Seat Height	440mm - 510mm	
Backrest Height	490mm	
Maximum Speed	10kph	
Range	54km	
Weight	110kg	
SWL	136kg	