



INVACARE PEGASUS PRO SCOOTER



SKU: 1637601-JR-H

PRODUCT DESCRIPTION

Purpose & Benefits of Invacare Pegasus Pro Scooter

The **Invacare Pegasus Pro Scooter** is a powerful, full-sized mobility scooter built for users who demand performance, comfort, and advanced safety features for long-distance travel. Designed for outdoor use across a range of terrains, the Pegasus Pro offers a smooth and stable ride with large pneumatic tyres, full suspension, and a high-capacity motor.

This scooter includes Invacare's intelligent electronics package, providing automatic speed reduction around corners and reliable diagnostics to enhance user confidence. With a high-back padded seat, adjustable tiller, and robust construction, the Pegasus Pro is ideal for active individuals who require reliable, long-range mobility support.

Key Features

- Full-sized scooter designed for extended outdoor use
- Powerful motor and battery for long-distance travel
- Full suspension system and large pneumatic tyres for superior ride comfort
- Automatic speed reduction for added safety while turning
- Ergonomic handlebars and adjustable tiller for personalised control

- High-back padded seat with swivel and armrests for comfort and access
- Front and rear lighting system for improved visibility
- Suitable for rougher terrains and longer journeys

Specifications

- Weight capacity: 160 kg
- Product weight: 136 kg
- Range: 52 km
- Backrest width: 510 mm
- Seat depth: 470 mm
- Seat to floor height: Min. 440 mm, Max. 510 mm
- Backrest height: Comfort: 500 mm, Deluxe: 490 mm, Premium: 630 mm
- Total length: 1320 mm
- Turning diameter: 2750 mm
- Battery: 2 x 75 Ah

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au