



## INVACARE PEGASUS PRO SCOOTER

**~~\$6,289.00~~ Original price was:**

**~~\$6,289.00.~~\$5,660.10 Current price is: \$5,660.10.**



**SKU: 1637601-JR-H**

## PRODUCT DESCRIPTION

### Purpose & Benefits of Invacare Pegasus Pro Scooter

The **Invacare Pegasus Pro Scooter** is a powerful, full-sized mobility scooter built for users who demand performance, comfort, and advanced safety features for long-distance travel. Designed for outdoor use across a range of terrains, the Pegasus Pro offers a smooth and stable ride with large pneumatic tyres, full suspension, and a high-capacity motor.

This scooter includes Invacare's intelligent electronics package, providing automatic speed reduction around corners and reliable diagnostics to enhance user confidence. With a high-back padded seat, adjustable tiller, and robust construction, the Pegasus Pro is ideal for active individuals who require reliable, long-range mobility support.

### Key Features

- Full-sized scooter designed for extended outdoor use
- Powerful motor and battery for long-distance travel
- Full suspension system and large pneumatic tyres for superior ride comfort
- Automatic speed reduction for added safety while turning
- Ergonomic handlebars and adjustable tiller for personalised control

- High-back padded seat with swivel and armrests for comfort and access
- Front and rear lighting system for improved visibility
- Suitable for rougher terrains and longer journeys

## Specifications

- Weight capacity: 160 kg
- Product weight: 136 kg
- Range: 52 km
- Backrest width: 510 mm
- Seat depth: 470 mm
- Seat to floor height: Min. 440 mm, Max. 510 mm
- Backrest height: Comfort: 500 mm, Deluxe: 490 mm, Premium: 630 mm
- Total length: 1320 mm
- Turning diameter: 2750 mm
- Battery: 2 x 75 Ah