




## INVACARE REA AZALEA TRANSIT WHEELCHAIR

**\$4,584.00**



**SKU:** AZALEA-TR-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Invacare Rea Azalea Transit Wheelchair is the natural choice for the moderately active user as well as for the more passive one. The comfortable seating system allows the user to alter the seat position, ensuring maximum comfort at all times. This combined with good driving characteristics makes Azalea the ideal choice. The new backrest, Flex 3, is multi-adjustable in height, width, and shape for different body forms. A self-propelled version is also available.

### FEATURES

- The Invacare Rea Azalea Transit Wheelchair features a Flex 3 tension adjustable backrest and backrest plates, providing superior customisation to ensure comfort and proper posture.
- The 25° seat tilt and 30° backrest recline allow users to find the most comfortable sitting position, reducing strain and optimising long-term usage.
- Unique to this model is the biomechanical positioning hinge on the backrest, promoting natural movement and increased mobility for the user.
- The wheelchair is available in 3 seating width options, offering a personalised fit to accommodate various body types.

### SPECIFICATIONS

- Width: seat width + 23cm
- Length: 102cm
- Product Weight: 34kg
- Warranty: 2 years
- Max User Weight: 135kg
- Seat Width: 40.6cm (16"); 45.72cm (18"); 50.8cm (20")
- Seat Depth: 43cm - 50cm

- Seat Height: 45cm
- Tilt Range: 25°

## **COMPLIANCE**

ARTG 99518

AS|NZS 3695

Crash Tested: ISO 7176-19

## **WARRANTY**

Frame: 2 years | Components: 2 years