



INVACARE REA AZALEA TRANSIT WHEELCHAIR



SKU: AZALEA-TR-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Invacare Rea Azalea Transit Wheelchair is the natural choice for the moderately active user as well as for the more passive one. The comfortable seating system allows the user to alter the seat position, ensuring maximum comfort at all times. This combined with good driving characteristics makes Azalea the ideal choice. The new backrest, Flex 3, is multi-adjustable in height, width, and shape for different body forms. A self-propelled version is also available.

FEATURES

- The Invacare Rea Azalea Transit Wheelchair features a Flex 3 tension adjustable backrest and backrest plates, providing superior customisation to ensure comfort and proper posture.
- The 25° seat tilt and 30° backrest recline allow users to find the most comfortable sitting position, reducing strain and optimising long-term usage.
- Unique to this model is the biomechanical positioning hinge on the backrest, promoting natural movement and increased mobility for the user.
- The wheelchair is available in 3 seating width options, offering a personalised fit to accommodate various body types.

SPECIFICATIONS

• Width: seat width + 23cm

• Length: 102cm

Product Weight: 34kgWarranty: 2 years

Max User Weight: 135kg

• Seat Width: 40.6cm (16"); 45.72cm (18"); 50.8cm (20")

• Seat Depth: 43cm - 50cm

• Seat Height: 45cm



• Tilt Range: 25°

COMPLIANCE

ARTG 99518 AS|NZS 3695

Crash Tested: ISO 7176-19

WARRANTY

Frame: 2 years | Components: 2 years

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 319903 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au