



INVACARE TOPAZ 9000 SELF PROPELLING WHEELCHAIR



SKU: 9TPZ

PRODUCT DESCRIPTION

The Invacare wheelchair has a variety of seat widths and depths, frame reinforcements that include gusseted sides, caster journals and dual heavy-duty cross braces, this chair provides not only ease in manoeuvrability but proven stability.

Key Features

- Available in vinyl or nylon upholstery
- Durable low-maintenance, carbon steel frame is long-lasting
- Reinforced frame includes gusseted sides, caster journals and dual heavy-duty cross braces
- Heavy-duty front riggings provide proper positioning for feet
- Heavy-duty arms designed to accommodate weight capacity and provide proper positioning
- Dual axle positions allow for quick and simple seat-to- oor height variations, 17.5" and 19.5"
- Heavy-duty package option for 1,000 lb. weight capacity

Seat width

• 20", 22", 24", 26", 28", 30"

Max. user weight

• 318kg 454kg (HD)

Seat depth

• 18", 20"

Seat height



• 15.5", 17.5", 19.5"

Backrest height

• 16", 17", 18"

Buy the Invacare heavy duty wheelchair online at MediMart using Afterpay or NDIS method today!

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au