



## INVACARE TOPAZ 9000 SELF PROPELLING WHEELCHAIR

**\$4,471.00** Original price was:

**~~\$4,471.00~~ \$4,023.90** Current price is: **\$4,023.90.**



**SKU: 9TPZ**

## PRODUCT DESCRIPTION

The Invacare wheelchair has a variety of seat widths and depths, frame reinforcements that include gusseted sides, caster journals and dual heavy-duty cross braces, this chair provides not only ease in manoeuvrability but proven stability.

## Key Features

- Available in vinyl or nylon upholstery
- Durable low-maintenance, carbon steel frame is long-lasting
- Reinforced frame includes gusseted sides, caster journals and dual heavy-duty cross braces
- Heavy-duty front riggings provide proper positioning for feet
- Heavy-duty arms designed to accommodate weight capacity and provide proper positioning
- Dual axle positions allow for quick and simple seat-to-floor height variations, 17.5" and 19.5"
- Heavy-duty package option for 1,000 lb. weight capacity

### Seat width

- 20", 22", 24", 26", 28", 30"

### Max. user weight

- 318kg 454kg (HD)

### Seat depth

- 18", 20"

### Seat height

- 15.5", 17.5", 19.5"

**Backrest height**

- 16", 17", 18"

Buy the Invacare heavy duty wheelchair online at MediMart using Afterpay or [NDIS method](#) today!