

JAY FUSION WHEELCHAIR CUSHION WITH FLUID INSERT

\$1,243.00 - \$1,597.00 Price range: \$1,243.00 through \$1,597.00



SKU: JFUSION1414-H



PRODUCT DESCRIPTION

Purpose & Benefits

The **Jay Fusion Wheelchair Cushion with Fluid Insert** delivers superior skin protection and stability for individuals at high risk of pressure injuries. Designed for extended wheelchair use, this cushion merges the skin-protecting features of the J2 with the exceptional postural support of the J3.

Its contoured closed-cell foam base, layered construction, and increased fluid volume promote optimal immersion, pressure distribution, and stability. The JAY Flow fluid tripad adjusts to the user's shape and sitting posture while maintaining proper positioning and skin integrity.

This cushion is ideal for users who experience weight changes, have limited mobility, or require consistent protection throughout the day. Whether for everyday use or high-dependency needs, the Fusion offers customisable support with excellent long-term performance.

Features

- Pre-contoured closed-cell foam base for enhanced stability
- JAY Flow fluid tripad insert for superior skin protection
- Adjustable positioning and pressure redistribution
- Suitable for individuals unable to weight shift independently

- Dual-cover system: stretch outer cover and water-resistant inner liner
- Aqua-guard zipper with anti-wicking thread
- 3DX spacer fabric promotes airflow and temperature regulation
- X-static® silver-infused fibres for antimicrobial protection
- Available with incontinence cover and AirExchange cover options
- Designed for extremely high-risk users and extended wheelchair use

Specifications

- **Insert Type:** JAY Flow™ Fluid Insert
- **Cushion Type:** Skin protection & positioning
- **Base:** Pre-contoured, closed-cell foam
- **Cover System:** Dual-layer, water-resistant inner + stretch outer
- **User Suitability:**
 - High risk of skin breakdown
 - Limited mobility or weight-shifting ability
 - Extended seating durations
 - Postural asymmetries
 - Weight fluctuation
- **Manufacturer:** Sunrise Medical
- **Model:** Jay Fusion

[Download The Jay Cushion Brochure Here](#)