

JAY J2 OBLIQUITY PAD

\$38.56



SKU: JBALCOP-ALH-H



PRODUCT DESCRIPTION

Purpose & Benefits

The **Jay J2 Obliquity Pad** is designed to enhance pelvic alignment and improve postural symmetry in wheelchair users who present with a fixed or flexible pelvic obliquity. When used with JAY cushions, particularly the Jay J2, it allows clinicians or users to accommodate or correct asymmetrical pelvic positioning for optimal stability, skin protection, and comfort.

The obliquity pad helps reduce the risk of pressure injuries caused by uneven weight distribution and is an essential accessory for fine-tuning seating systems for users with complex positioning needs.

Features

- Customisable wedge pad for managing pelvic obliquity
- Designed to be placed under the cushion's pelvic loading area
- Can be used to accommodate or correct asymmetry
- Compatible with **Jay J2** and other JAY cushions
- Provides enhanced pelvic alignment and skin protection
- Lightweight, low-profile, and easy to position

- Ideal for wheelchair users with postural asymmetries

Specifications

- **Product Type:** Seating accessory
- **Compatibility:** Jay J2 and most JAY cushions
- **Application:** Pelvic obliquity accommodation or correction
- **Placement:** Under cushion base or internal layer depending on setup
- **Material:** High-density foam
- **Manufacturer:** Sunrise Medical
- **Model:** J2 Obliquity Pad