


KARMA AGILE SELF PROPELLED WHEELCHAIR



SKU: KMAGL-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

It's all about performance, comfort, and style. Meet Agile (AGL), the new standard for transferring wheelchairs.

Agile not only inherited its multi-adjustability, but also evolved to provide a wider variety of options to serve each individual's needs.

The Karma Agile economy wheelchair has a sleek design that comes in several sizes and is packed with high quality features you expect from Karma. The Agile's overall weight has been reduced for easy lifting and storage.

FEATURES

- Lightweight design for easy manoeuvrability
- Foldable frame for convenient storage and transportation
- Ergonomic seating for comfort and support
- Quick-release rear wheels for easier transport
- Adjustable footrests for customised comfort
- Puncture-proof tires for worry-free use
- Available in a range of sizes to suit different user needs
- Stylish design with a variety of colour options
- Designed for both indoor and outdoor use.

SPECIFICATIONS

- **Seat Width Options:**

- 40.6cm (16")
- 45.72cm (18")
- 50.8cm (20")
- **Seat Depth Options:**
 - 40.6cm (16")
 - 45.72cm (18")
- **Frame Colour Options:**
 - Black
 - Silver
- **Fabric Colour Options:**
 - Black
 - Grey Dot
- **Seat Height:** 50cm
- **Back Height:** 43cm
- **Safe Working Load:** 125kg
- **Weight:** 14kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au