



KARMA ERGO LITE 2 DELUXE SELF PROPELLED WHEELCHAIR



SKU: KM2512F-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Karma Ergo Lite 2 Deluxe Self Propelled Wheelchair is a lightweight and durable mobility aid designed to enhance user comfort and convenience. With a weight of only 10.3 kg, this wheelchair is easy to handle and transport, making it an excellent choice for active users who are always on the go.

The Ergo Lite 2 Deluxe features a unique S-shaped ergonomic seat that promotes a better seating posture and reduces pressure on the lower back. The backrest and footrest can be adjusted to fit the user's preferences, ensuring a comfortable and customised fit.

FEATURES

- Lightweight Aluminium Folding Frame
- Swing-away and Removable Footrests
- Ergonomic Pad Armrest
- Fixed 20" x 1 3/8 Solid Tyre Wheels
- 6" Solid Castors
- Push to Lock Brakes
- Fixed Seat Height
- Folding Back
- Dacron Fire Retardant Upholstery
- 2-Year Warranty on Frame

SPECIFICATIONS

- **Seat Width:** 40.6cm (16"); 45.72cm (18");
- **Seat Depth:** 40.6cm (16")
- **Total Width:** 65cm
- **Total Length:** 93cm
- **Total Height:** 91cm
- **Width (folded):** 31cm
- **Armrest Height:** 21.5cm
- **Seat Height:** 48cm
- **Backrest Height:** 42cm
- **Product Weight:** 10.3kg
- **Maximum User Weight:** 100kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au

