



KARMA ERGO LITE 2 DELUXE TRANSIT WHEELCHAIR

\$1,325.00



SKU: 2-KM2512F1-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Karma Ergo Lite 2 Deluxe Transit Wheelchair is a lightweight and versatile mobility aid that is perfect for individuals who require assistance when moving around. This transit wheelchair features a comfortable and ergonomic seat that ensures proper posture and reduces pressure on the lower back, making it ideal for long-term use.

The wheelchair is equipped with large rear wheels, making it easy to manoeuvre and providing a smooth and comfortable ride. The deluxe version of the Ergo Lite 2 Transit Wheelchair also features an adjustable backrest, allowing users to find the perfect level of support for their back.

In addition to its ergonomic design, the Ergo Lite 2 Deluxe Transit Wheelchair is also easy to transport and store, thanks to its lightweight construction and folding frame.

FEATURES

- Lightweight Aluminium Folding Frame
- Swing-away and Removable Footrests
- Fixed 12" x 2" Solid Tyre Wheels
- 6" Solid Castors
- Attendant Brakes
- Fixed Seat Height
- Folding Back

SPECIFICATIONS

Sizes (width x depth)	41cm x 41 cm (16"x16")	46cm x 41cm (18"x16")
Height	91cm	91cm
Weight	8.7kg	8.8kg
Total Width	63cm	65cm
Total Length	87cm	87cm
Seat Height	48cm	48cm
Backrest Height	42cm	42cm
Width (Folded)	31cm	31cm
Maximum User Weight	100kg	100kg