




KARMA ERGO LITE 2 DELUXE TRANSIT WHEELCHAIR



SKU: KM2512F1-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Karma Ergo Lite 2 Deluxe Transit Wheelchair is a lightweight and versatile mobility aid that is perfect for individuals who require assistance when moving around. This transit wheelchair features a comfortable and ergonomic seat that ensures proper posture and reduces pressure on the lower back, making it ideal for long-term use.

The wheelchair is equipped with large rear wheels, making it easy to manoeuvre and providing a smooth and comfortable ride. The deluxe version of the Ergo Lite 2 Transit Wheelchair also features an adjustable backrest, allowing users to find the perfect level of support for their back.

In addition to its ergonomic design, the Ergo Lite 2 Deluxe Transit Wheelchair is also easy to transport and store, thanks to its lightweight construction and folding frame.

FEATURES

- Lightweight Aluminium Folding Frame
- Swing-away and Removable Footrests
- Fixed 12" x 2" Solid Tyre Wheels
- 6" Solid Castors
- Attendant Brakes
- Fixed Seat Height
- Folding Back

SPECIFICATIONS

Sizes (width x depth)	41cm x 41 cm (16"x16")	46cm x 41cm (18"x16")
Height	91cm	91cm
Weight	8.7kg	8.8kg
Total Width	63cm	65cm
Total Length	87cm	87cm
Seat Height	48cm	48cm
Backrest Height	42cm	42cm
Width (Folded)	31cm	31cm
Maximum User Weight	100kg	100kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
 1300 633 468

support@medimarnt.com.au