



## **KARMA ERGO LITE 2 TALL SELF-PROPEL WHEELCHAIR 20"**



**SKU:** KM2512TLF22-20

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits**

The **Karma Ergo Lite 2 Tall Self-Propel Wheelchair - 20"** is a lightweight, foldable wheelchair designed for users who prefer independent mobility. Featuring Karma's patented S-Ergo seating system, it promotes comfort and proper posture while remaining easy to transport and store. With an increased seat height and width, it accommodates taller users and provides smooth manoeuvrability with large self-propelling rear wheels.

### **Features**

- Lightweight folding frame for effortless transport and storage
- Self-propel design for independent use and smooth operation
- S-Ergo patented seating system encourages correct posture and long-term comfort
- Push and pull attendant brake system for added safety and control
- Detachable swing-away footrests for convenient transfers
- Armrests designed to allow closer positioning to desks or tables
- Top backrest section folds down for compact transportation

- Crash tested and compliant with ISO 7176-19 for transport safety

## Specifications

Specification	Detail
Overall Length	1030 mm
Overall Width	700 mm
Overall Height	950 mm
Seat Width	510 mm (20")
Seat Depth	435 mm
Floor to Seat Height	510 mm
Folded Width	280 mm
Folded Depth	755 mm
Armrest Height	210 mm
Backrest Height	435 mm
Product Weight	11.8 kg
Maximum User Weight	115 kg
Brand	K Care

**Note:** Also available in a transit (attendant-propelled) model. Please select the correct configuration to suit your mobility needs.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

**Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

**Monash Showroom**

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

**NT**

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)