



KARMA ERGO LITE 2 TALL TRANSIT WHEELCHAIR 20"



SKU: KM2512TLF14-20

PRODUCT DESCRIPTION

Purpose & Benefits

The **Karma Ergo Lite 2 Tall Transit Wheelchair - 20"** is a lightweight and foldable attendant-propelled wheelchair designed for comfort, practicality, and ease of use. Built with the patented S-Ergo seating system, it promotes proper posture and long-term comfort. Its increased seat height, width, and weight capacity accommodate a wider range of users, making it an ideal choice for daily mobility and transport.

Features

- Lightweight folding frame for easy handling and storage
- Patented S-Ergo seating system promotes correct posture and comfort
- Increased seat width and floor-to-seat height for taller users
- Push and pull attendant brake system for added control and safety
- Detachable swing-away footrests for convenient transfers
- Armrest design allows closer positioning to desks or tables
- Top backrest section folds down for compact transportation

- Crash tested in accordance with ISO 7176-19 standards

Specifications

Specification	Detail
Overall Length	1030 mm
Overall Width	665 mm
Overall Height	950 mm
Seat Width	510 mm (20")
Seat Depth	435 mm
Floor to Seat Height	510 mm
Folded Width	270 mm
Folded Depth	720 mm
Armrest Height	210 mm
Backrest Height	435 mm
Product Weight	10.3 kg
Maximum User Weight	115 kg
Brand	Karma

Note: Available in transit (attendant-propelled) and self-propelled models. Please confirm your preferred configuration before ordering.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au