



KARMA ERGO LITE TRANSIT WHEELCHAIR



SKU: KM2501F-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Karma Ergo Lite Transit Wheelchair is a lightweight and compact transit wheelchair designed for easy transportation and storage. It features an ergonomic S-shaped seat that conforms to the natural shape of the body, providing maximum comfort and support. The chair is made of high-grade aluminium, which makes it incredibly lightweight yet strong and durable. The backrest is adjustable and can be tilted to fit the user's preference, and the footrests can be adjusted in height to provide a comfortable and supportive footrest. The Ergo Lite Transit Wheelchair also features easy-to-use brakes, making it simple for the user or caregiver to control the chair's speed and direction. With a weight capacity of up to 100 kg, this wheelchair is perfect for those who need a lightweight and easily transportable chair for everyday use.

FEATURES

- **Ultra-Lightweight Design:** Weighing only 8.3 kg, this wheelchair is one of the lightest transit chairs available, making it easy to handle and transport.
- **Ergonomic Seating System:** Features the patented S-Ergo seating system, which is designed to reduce pressure on the user's body, enhance comfort, and prevent sliding.
- **Height Adjustable Footrest:** The footrests are adjustable to suit various user heights, ensuring maximum comfort and support during transit.
- **Fixed Ergonomic Armrest Pad:** Includes fixed armrests with ergonomic padding, providing comfortable support for the user's arms during use.
- **Attendant Braking System:** The wheelchair is equipped with attendant-operated brakes, allowing the caregiver to safely control the chair during movement.
- **Dacron, Fire-Retardant Upholstery:** Upholstered with high-quality Dacron fabric that is fire

retardant, adding an extra layer of safety for users.

- **Foldable for Easy Storage:** The Karma Ergo Lite can be quickly folded, making it easy to store in small spaces or transport in vehicles.
- **Puncture-Proof Tyres:** Equipped with solid, puncture-proof tyres for worry-free use on various terrains, offering a smooth ride indoors and outdoors.

SPECIFICATIONS

Sizes	41x41cm (16x16")	46x41cm (18x16")
Total Height	91cm	91cm
Total Length	94cm	94cm
Total Width	57cm	62cm
Seat Width	40cm	45cm
Seat Height	48cm	48cm
Armrest Height	21.5cm	21.5cm
Backrest Height	42cm	42cm
Seat Depth	42cm	42cm
Width (Folded)	27.5cm	27.5cm
Product Weight	8.1kg	8.3kg
Maximum User Weight	100kg	100kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au