




KARMA FLEXX HD SELF PROPELLED WHEELCHAIR



SKU: N/A

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Flexx Series is all about comfort, support, and flexibility. The Flexx's elliptical-frame design is durable, stylish, and robust, supporting a wide array of adjustable features, such as rear-wheel axle position, seat depth, backrest angle, and backrest height. The Flexx HD is a heavy-duty version with larger seat options and a safe working load of 170kg.

KEY FEATURES

- Height adjustable flip back padded armrests
- Detachable swing away height adjustable footrests
- Multiple self-propel wheel settings
- Patented ergonomic hand rim provide extra grip and easy propulsion
- Optional stump supports and elevating leg rests
- Crash test approved in accordance with ISO 7176-19 only in 16" and 18" seat width models. The Flexx HD models are not crash test approved.

SPECIFICATIONS

Feature	20" Wheelchair	22" Wheelchair	24" Wheelchair
Overall Length	1180 mm	1180 mm	1180 mm
Overall Width	810 mm	810 mm	810 mm
Overall Height	1040 mm	1040 mm	1040 mm
Seat Width	510 mm	560 mm	610 mm
Seat Depth	410 mm	460 mm	510 mm
Floor to Seat Height	520 mm	520 mm	520 mm

Feature	20" Wheelchair	22" Wheelchair	24" Wheelchair
Folded Width	350 mm	350 mm	350 mm
Folded Depth	1180 mm	1180 mm	1180 mm
Folded Height	1040 mm	1040 mm	1040 mm
Min Armrest Height	220 mm	220 mm	220 mm
Max Armrest Height	300 mm	300 mm	300 mm
Backrest Height	430 mm	430 mm	430 mm
Product Weight	16.5 kg	17 kg	17.6 kg
Max User Weight	170 kg	170 kg	170 kg
Brand	Karma	Karma	Karma

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
 1300 633 468

support@medimarnt.com.au