



KARMA FLEXX SELF PROPELLED WHEELCHAIR

\$1,829.00 - \$2,587.00Price range: \$1,829.00 through \$2,587.00



SKU: KM8022-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Karma Flexx Self Propelled Wheelchair is a versatile and lightweight mobility aid designed for people with mobility impairments. Its unique design allows for easy customisation and adjustment, making it comfortable and suitable for a wide range of users. The chair is easy to manoeuvre and provides excellent stability, making it ideal for indoor and outdoor use. The Flexx Self Propelled Wheelchair is equipped with a number of features to enhance the user's experience, including adjustable armrests and footrests, a comfortable seat cushion, and a foldable frame for easy storage and transportation.

FEATURES

- Folding Lightweight Aluminium Frame
- Removable, Height Adjustable, Swing in and out Footrest
- Angle Adjustable Backrest
- Swing away, Height Adjustable with Ergonomic Pad Armrests
- Quick release Wheels
- 7" x 1" Castors
- Push to Lock Brakes
- Depth Height Adjustable Seat
- Dacron, Fire Retardant Upholstery
- 5-Year Warranty

SPECIFICATIONS

· Seat width



- 40.6cm (16")
- 45.72cm (18")
- o 50.8cm (20")
- o 56cm (22")
- Seat depth
 - 40.6cm (16")
 - 45.72cm (18")
 - ∘ 50.8cm (20")
- Seat height (from floor) 390-530mm
- Seat height (from footrest) 380-480mm
- Backrest height 380-530mm
- Backrest adjustment angle -4/+16 degrees
- Folded width 340mm
- Overall weight
 - 18" option 15.9kg
 - ∘ 22" option -18kg
- Wheel size
 - ∘ Front 7"
 - Rear 24"
- Maximum weight capacity 130kg