


KARMA FLEXX SELF PROPELLED WHEELCHAIR



SKU: KM8022-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Karma Flexx Self Propelled Wheelchair is a versatile and lightweight mobility aid designed for people with mobility impairments. Its unique design allows for easy customisation and adjustment, making it comfortable and suitable for a wide range of users. The chair is easy to manoeuvre and provides excellent stability, making it ideal for indoor and outdoor use. The Flexx Self Propelled Wheelchair is equipped with a number of features to enhance the user's experience, including adjustable armrests and footrests, a comfortable seat cushion, and a foldable frame for easy storage and transportation.

FEATURES

- Folding Lightweight Aluminium Frame
- Removable, Height Adjustable, Swing in and out Footrest
- Angle Adjustable Backrest
- Swing away, Height Adjustable with Ergonomic Pad Armrests
- Quick release Wheels
- 7" x 1" Castors
- Push to Lock Brakes
- Depth Height Adjustable Seat
- Dacron, Fire Retardant Upholstery
- 5-Year Warranty

SPECIFICATIONS

- Seat width

- 40.6cm (16")
- 45.72cm (18")
- 50.8cm (20")
- 56cm (22")
- Seat depth
 - 40.6cm (16")
 - 45.72cm (18")
 - 50.8cm (20")
- Seat height (from floor) – 390-530mm
- Seat height (from footrest) – 380-480mm
- Backrest height – 380-530mm
- Backrest adjustment angle – -4/+16 degrees
- Folded width – 340mm
- Overall weight
 - 18" option – 15.9kg
 - 22" option – 18kg
- Wheel size
 - Front 7"
 - Rear 24"
- Maximum weight capacity – 130kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au