



KARMA S-ERGO 125 SELF PROPELLED WHEELCHAIR

\$1,606.00 - \$1,942.00 Price range: \$1,606.00 through \$1,942.00



SKU: KMS125

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Karma S-Ergo 125 Self Propelled Wheelchair is a top-of-the-line mobility aid designed for comfort, style, and functionality. With its ergonomic design and lightweight frame, this wheelchair is perfect for active individuals who need a reliable and comfortable way to get around. Its advanced features include adjustable seat height, quick-release wheels, and a foldable design for easy storage and transport. Whether you're using it for daily activities or for travel, the Karma S-Ergo 125 Self Propelling Wheelchair is an excellent choice for anyone looking for a high-quality wheelchair that combines form and function.

FEATURES

- Lightweight Aluminium Frame
- Swing-away Removable Footrest
- Swing-away Armrest with Ergonomic Pads
- Quick Release Solid Tyre Wheels
- Solid 7" Castors
- Push-to-lock with Attendant Braking System
- Fixed Seat Height
- Folding Back
- Mesh with Aegis Microbe Shield Padding Upholstery

SPECIFICATIONS

Seat Width	41cm (16")	46cm (18")	51cm (20")
------------	------------	------------	------------

Total Width	59cm	64cm	69cm
Total Length	100cm	100cm	100cm
Total Height	91cm	91cm	91cm
Width (Folded)	30cm	30cm	30cm
Seat Height	51cm	51cm	51cm
Armrest Height	20cm	20cm	20cm
Backrest Height	40cm	40cm	40cm
Product Weight	14.2kg	14.5kg	14.8kg
Maximum User Weight	115kg	115kg	115kg