




KARMA S-ERGO 125 TRANSIT WHEELCHAIR

\$1,445.40 - \$1,747.80 Price range: \$1,445.40 through \$1,747.80



SKU: 2-KMS125T-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Karma S-Ergo 125 Transit is equipped with the Unique S-Ergo (Smart Ergo) S shaped seat.

The S-Ergo design was developed over two years, after extensive testing by Karma's Research and Development Department. The result is a leap forward in wheelchair design.

In Addition to the S-Ergo seating system, the S-Ergo 125 Transit also has seat upholstery and pads treated with Aegis Microbe shield, to reduce stains and odour. The streamline oval tubes not only make the wheelchair look great, but also add to the strength of the wheelchair structure, which is incredibly light. Frame made from aircraft aluminium.

The S-Ergo 125 Transit also features patented footrests, which swing inward and outward for added convenience.

FEATURES

- Folding Lightweight Aluminium Frame
- Swing-away Removable Footrest
- Swing-away with Ergonomic Pads Armrest
- Fixed Solid Tyre Wheels 14" x 1.75
- Solid 7" Castors
- Push-to-lock with Attendant Braking System
- Fixed Seat Height
- Folding Back
- Mesh with Aegis Microbe Shield Padding Upholstery

SPECIFICATIONS

Sizes	41x43cm (16x17")	46x43cm (18x17")	51x43cm (20x17")
Total Width	63cm	68cm	73cm
Total Length	108cm	108cm	108cm
Total Height	97cm	97cm	97cm
Armrest Height	20cm	20cm	20cm
Seat Height	50cm	50cm	50cm
Seat Width	40cm	45cm	50cm
Seat Depth	43cm	43cm	43cm
Backrest Height	40cm	40cm	40cm
Width (Folded)	31cm	31cm	31cm
Product Weight	13.1kg	13.1kg	13.3kg
Maximum User Weight	115kg	115kg	115kg