



KARMA TRANSIT II WITH FOLDING BACK TRANSIT WHEELCHAIR



SKU: KM8520F-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Karma Transit II with Folding Back Transit Wheelchair is a reliable and practical mobility aid that is ideal for users who require assistance while travelling. This lightweight wheelchair is equipped with a foldable backrest, making it easy to store and transport. The Transit II features a comfortable seat and a padded backrest, as well as detachable and adjustable footrests to ensure maximum comfort for the user.

The Transit II also comes with a range of safety features, including locking brakes and anti-tipping devices, which provide peace of mind and security for the user. The wheelchair is constructed with durable materials, ensuring that it can withstand daily use and last for years to come.

FEATURES

- A lightweight, durable, and reliable manual wheelchair that makes it the right choice for both user and carer
- Folding, lightweight aluminium frame with foldable backrest for ease of transfers
- Ideal for those who travel with assistance
- Offers an impressive range of standard features with an extensive selection of personal customising available
- Aircraft grade materials and unique designs of crossbars and link bars provide strength
- Detachable and height adjustable armrests and footrests, plus pull to lock brakes

SPECIFICATIONS

| Spec | 40.6 cm x 40.6 cm (16"x16") | 45.7 cm x 40.6 cm (18" x 16") | 45.7 cm x 45.7 cm (18" x 18") | 50.8 cm x 45.7 cm (20" x 18") | 55.9 cm x 45.7 cm (22" x 18") |
|----------------------------|-----------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| Seat Width | 41cm | 46cm | 46cm | 51cm | 55cm |
| Seat Depth | 41cm | 41cm | 46cm | 46cm | 46cm |
| Total Width | 60cm | 65cm | 65cm | 69cm | 74cm |
| Total Length | 105cm | 105cm | 105cm | 105cm | 105cm |
| Total Height | 92cm | 92cm | 92cm | 92cm | 92cm |
| Armrest Height | 23-28cm | 23-28cm | 23-28cm | 23-28cm | 23-28cm |
| Seat Height | 52cm | 52cm | 52cm | 52cm | 52cm |
| Backrest Height | 41-47cm | 41-47cm | 41-47cm | 41-47cm | 41-47cm |
| Product Weight | 13.5kg | 14.3kg | 14.3kg | 15.2kg | 15.7kg |
| Maximum User Weight | 130kg | 130kg | 130kg | 160kg | 160kg |

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au