



## **KARMA TRANSIT II WITH FOLDING BACK TRANSIT WHEELCHAIR**

**\$1,718.00 - \$2,797.00** Price range: \$1,718.00 through \$2,797.00



**SKU:** KM8520F-H

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

The Karma Transit II with Folding Back Transit Wheelchair is a reliable and practical mobility aid that is ideal for users who require assistance while travelling. This lightweight wheelchair is equipped with a foldable backrest, making it easy to store and transport. The Transit II features a comfortable seat and a padded backrest, as well as detachable and adjustable footrests to ensure maximum comfort for the user.

The Transit II also comes with a range of safety features, including locking brakes and anti-tipping devices, which provide peace of mind and security for the user. The wheelchair is constructed with durable materials, ensuring that it can withstand daily use and last for years to come.

### **FEATURES**

- A lightweight, durable, and reliable manual wheelchair that makes it the right choice for both user and carer
- Folding, lightweight aluminium frame with foldable backrest for ease of transfers
- Ideal for those who travel with assistance
- Offers an impressive range of standard features with an extensive selection of personal customising available
- Aircraft grade materials and unique designs of crossbars and link bars provide strength
- Detachable and height adjustable armrests and footrests, plus pull to lock brakes

## **SPECIFICATIONS**

<b>Spec</b>	<b>40.6 cm x 40.6 cm (16"x16")</b>	<b>45.7 cm x 40.6 cm (18" x 16")</b>	<b>45.7 cm x 45.7 cm (18" x 18")</b>	<b>50.8 cm x 45.7 cm (20" x 18")</b>	<b>55.9 cm x 45.7 cm (22" x 18")</b>
<b>Seat Width</b>	41cm	46cm	46cm	51cm	55cm
<b>Seat Depth</b>	41cm	41cm	46cm	46cm	46cm
<b>Total Width</b>	60cm	65cm	65cm	69cm	74cm
<b>Total Length</b>	105cm	105cm	105cm	105cm	105cm
<b>Total Height</b>	92cm	92cm	92cm	92cm	92cm
<b>Armrest Height</b>	23-28cm	23-28cm	23-28cm	23-28cm	23-28cm
<b>Seat Height</b>	52cm	52cm	52cm	52cm	52cm
<b>Backrest Height</b>	41-47cm	41-47cm	41-47cm	41-47cm	41-47cm
<b>Product Weight</b>	13.5kg	14.3kg	14.3kg	15.2kg	15.7kg
<b>Maximum User Weight</b>	130kg	130kg	130kg	160kg	160kg