

## **KARMA VIP 2 TILT/RECLINE SELF-PROPELLED WHEELCHAIR**



**SKU:** KMVIP2-H



### **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

The Karma VIP 2 Tilt/Recline Self-Propelled Wheelchair is a lightweight, foldable, and durable wheelchair that offers the best of both worlds: a tilt-in-space mechanism for pressure relief and a reclining backrest for comfort.

The tilt-in-space mechanism allows the entire seat and backrest to tilt backward up to 35 degrees, which helps to relieve pressure on the hips and buttocks. This can help to prevent the development of pressure sores and improve circulation. The reclining backrest can be adjusted to a variety of positions, from upright to fully reclined, so users can find the most comfortable position for them.

The Karma VIP 2 is also a self-propelled wheelchair, so users can maintain their independence and mobility. It has a lightweight frame that makes it easy to push and manoeuvre, and it folds up for easy transport.

### **FEATURES**

- Tilt-in-space and anti-sliding reclining system features anti-shear.
- Allows users with partial trunk control to traverse easily and in complete comfort.
- Lightweight, foldable and durable design.
- Height adjustable and detachable padded armrests.
- Detachable swing away footrests.
- Curved multi-adjustable padded headrest.

- AEGIS microbe shield inhibits bacterial and fungal growths, eliminates odours and is easy to clean.
- Tension adjustable backrest.
- Colour-coded features for intuitive adjustment.
- Modularised design allows easy upgrades, such as Jay and Matrix backrest systems.
- Crash test approved in accordance with ISO 7176-19.
- Wide range of options and accessories available.

The Karma VIP 2 Tilt/Recline Self-Propelled Wheelchair is a versatile and comfortable wheelchair that is ideal for users who need a combination of pressure relief and comfort. It is lightweight, foldable, and easy to use, making it a great choice for people who want to maintain their independence and mobility.

Buy this Karma VIP 2 Tilt/Recline Self-Propelled Wheelchair online with Medimart using Afterpay or [NDIS method](#) today!

## SPECIFICATIONS

Size	41cm (16")	46cm (18")
<b>Seat Width</b>	41cm	46cm
<b>Seat Depth</b>	46cm	46cm
<b>Total Height</b>	98cm	98cm
<b>Total Length</b>	105cm - 110cm	105cm - 110cm
<b>Total Width</b>	63.5cm - 69cm	63.5cm - 69cm
<b>Tilt Range</b>	0-35	0-35
<b>Recline Range</b>	0-30	0-30
<b>Floor to Seat Height</b>	48cm - 50cm	48cm - 50cm
<b>Product Weight</b>	23.4kg	23.4kg
<b>Maximum User Weight</b>	115kg	115kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
 Brighton SA 5048  
 08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
 Ave, Elizabeth Vale SA 5112  
 08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriiffs Rd,  
 Reynella SA 5161  
 08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)