




**KARMA VIP2 TILT & RECLINE TRANSIT
WHEELCHAIR 18"X16" - ELEVATED LEG RESTS**

\$3,728.00



SKU: K MVP ZT-18

PRODUCT DESCRIPTION

Purpose & Benefits

The **Karma VIP2 Tilt & Recline Transit Wheelchair** in 18" × 16" seat size with elevated leg rests gives users fluid mobility and premium comfort. Its tilt-in-space mechanism and anti-sliding reclining system minimise skin shear and maintain posture for users with partial trunk control. Lightweight yet durable and foldable, this wheelchair pairs positioning support with practicality—making transfers easier, ensuring pressure relief, and enhancing overall comfort during use.

Features

- Tilt-in-space and anti-sliding reclining system to reduce shear and maintain spinal alignment
- Lightweight, foldable aluminium frame for easy transport and storage
- Height-adjustable, fully detachable, padded armrests—no tools needed
- Detachable swing-away footrests for expedient transfers
- Curved adjustable headrest and tension-adjustable backrest for enhanced positioning
- Crash tested to ISO 7176-19 for safe occupant transportation

Specifications

- **Overall Length:** 1130 mm
- **Overall Width:** 680 mm
- **Overall Height:** 1100 mm
- **Seat Width:** 450 mm (18")
- **Seat Depth:** 415 mm (16")
- **Floor-to-Seat Height:** 490 mm
- **Armrest Height:** 270 mm
- **Backrest Height:** 490 mm
- **Folded Dimensions (W × D × H):** 420 × 840 × 700 mm
- **Product Weight:** 17.7 kg
- **Maximum User Weight:** 115 kg