




REDGUM KATIE HIGH BACK ADJUSTABLE UTILITY CHAIR



SKU: N/A

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Katie High Back Chair is a modern adjustable chair featuring a contoured backrest for optimal postural support. A fantastic day chair for use at home and in-care facilities.

FEATURES

- Modern design styling Adjustable seat height
- Adjustable seat depth
- Contoured padded backrest
- Comfortable padded seat
- Large padded armrests
- Chair mobility wheels on rear legs as standard Backrest folds down for transit + storage
- Strong lightweight construction
- Epoxy powder coated framework
- Vinyl covered upholstery (various colours)

SPECIFICATIONS

Overall Width:	66.5cm
Overall Depth:	65cm
Overall Height:	114cm
Armrest Height:	22cm
Backrest Width:	48cm
Backrest Height:	70cm

Seat Width (between Armrests):	48cm
Seat Depth:	43.5 to 47.5cms
Seat Height:	45 to 57.5cms
Backrest Angle:	98°
Adjustment Increments:	2.5cm
Weight:	13kg
Weight Capacity:	200kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au