



REDGUM KATIE LOW BACK CHAIR

\$325.00



SKU: RG635

PRODUCT DESCRIPTION

This Katie low back chair supports the lower to the middle of the back. The backrest of the chair can be collapsed down for transportation + storage and is strong and lightweight.

Please note: This product is typically dispatched directly from the supplier. **Extended dispatch times** may apply. We'll contact you after purchase with an estimated delivery timeframe.

Key Features Include:

- Chair mobility wheels on rear legs as standard
- The backrest of the chair can be folded down for transit + storage
- Strong lightweight construction
- [Epoxy powder](#) coated framework
- Vinyl-covered upholstery (various colours)
- Adjustable seat depth
- Contoured padded backrest
- Comfortable padded seat
- Large padded armrests

This Katie low back chair supports the lower to the middle of the back. The backrest of the chair can be collapsed down for transportation+ storage and is strong and lightweight.

Key Features Include:

- Chair mobility wheels on rear legs as standard

- The backrest of the chair can be folded down for transit + storage
- Strong lightweight construction
- [Epoxy powder](#) coated framework
- Vinyl-covered upholstery (various colours)Adjustable seat depth
- Contoured padded backrest
- Comfortable padded seat
- Large padded armrests

SPECIFICATIONS

- Overall Dimensions (at lowest setting): 665mm (W) x 650mm (D) x 915mm (H)
- Seat Width: 480mm (between armrests)
- The Seat Depth is: 435 – 475mm
- Seat Height: 450 – 575mm (from floor to top of seat) 25mm adjustment increments
- Armrest Height: 220mm (from top of seat to top of armrest)
- The Backrest Width is : 480mm (at widest point)
- Backrest Height: 480mm
- Backrest Angle: 98°
- Mobility Wheel: Ø50mm
- Product Weight: 11kg
- Maximum User Weight: 200kg

Buy this Low Back Chair online at MediMart using Afterpay or [NDIS method](#) today!