



## **REDGUM KATIE LOW BACK CHAIR**

\$325.00 Original price was: \$325.00.\$293.00 Current price is: \$293.00.Sale Available Online Only!



**SKU: RG635** 

#### PRODUCT DESCRIPTION

This Katie low back chair supports the lower to the middle of the back. The backrest of the chair can be collapsed down for transportation + storage and is strong and lightweight.

**Please note:** This product is typically dispatched directly from the supplier. **Extended dispatch times** may apply. We'll contact you after purchase with an estimated delivery timeframe.

# **Key Features Include:**

- Chair mobility wheels on rear legs as standard
- The backrest of the chair can be folded down for transit + storage
- Strong lightweight construction
- **Epoxy powder** coated framework
- Vinyl-covered upholstery (various colours)
- Adjustable seat depth
- Contoured padded backrest
- · Comfortable padded seat
- Large padded armrests

This Katie low back chair supports the lower to the middle of the back. The backrest of the chair can be collapsed down for transportation+ storage and is strong and lightweight.

# **Key Features Include:**

· Chair mobility wheels on rear legs as standard



- The backrest of the chair can be folded down for transit + storage
- Strong lightweight construction
- **Epoxy powder** coated framework
- Vinyl-covered upholstery (various colours)Adjustable seat depth
- Contoured padded backrest
- Comfortable padded seat
- Large padded armrests

## **SPECIFICATIONS**

- Overall Dimensions (at lowest setting): 665mm (W) x 650mm (D) x 915mm (H)
- Seat Width: 480mm (between armrests)
- The Seat Depth is: 435 475mm
- Seat Height: 450 575mm (from floor to top of seat) 25mm adjustment increments
- Armrest Height: 220mm (from top of seat to top of armrest)
- The Backrest Width is: 480mm (at widest point)
- Backrest Height: 480mm
- Backrest Angle: 98°
- Mobility Wheel: Ø50mm
- Product Weight: 11kg
- Maximum User Weight: 200kg

Buy this Low Back Chair online at MediMart using Afterpay or NDIS method today!