




## KCARE ALIVIO JANNALI HOMECARE BED



**SKU:** KA3200KSBLK

## PRODUCT DESCRIPTION

### Purpose & Benefits of the Kcare Alivio Jannali Homecare Bed

Jannali's customisable positioning features enable you to find your optimal sleeping or resting position – with five (5) different pre-set positions and two memory functions using the Bluetooth phone app (iOS & Android compatible) to store your preferred bed positions.

An Autoregression section is built into the bed base to decrease pressure on the abdominal area and reduce shear forces on the occupant. The Jannali also features Trendelenburg and reverse Trendelenburg positioning.

The electric height adjustment function provides independence and reduces reliance on carers for assistance. For occupants requiring assisted transfers, it provides the carer with flexibility of height depending on the transfer equipment being used.

For additional safety, a non-entrapment zone under the bed stops the bed frame from touching the floor. The standard under-bed dimming LED light (*motion detection version optional*) provides visibility at night.

The full-featured Jannali has various convenient features, including a standard wired USB handset (or downloadable phone app), battery backup, and USB ports for charging devices.

Optional: The Jannali also comes with an exciting range of options, including matching headboard/footboard, side/corner mattress retainers, armrest/armrest cover, short/long bedrails, side pocket, height-adjustable bed stick, height-adjustable assist rail, and motion detection LED light.

Jannali comes in a range of modern, stylish colours and fabrics, including a soft-touch, easy-to-clean brushed vinyl (Black) and soft leather-look PU option (Beige and Canley), with colour availability for optional matching companion beds.

## Features

- A multi-function bed designed to enhance comfort and support.
- Electrically adjustable height, featuring one of the lowest platform heights available.
- Patented auto-regression backrest for optimal positioning.
- Trendelenburg and reverse Trendelenburg positioning (0-14°).
- Adjustable backrest angle (0-75°).
- Adjustable legrest angle (0-45°).
- Zero gravity pre-set function for ultimate relaxation.
- Additional memory setting to save user-specific positions.
- Home function lowers the bed to a flat position at its lowest height.
- Under-bed LED safety light, with optional motion detection.
- Standard backup battery for added peace of mind.
- A variety of rail and support stick options available.
- Optional head and footboard.
- Independently tested to comply with AS/NZS 60335.2.32 in conjunction with AS/NZS 60335.1.

## Colour Options

## Specifications

Specification	Long Single	King Single
<b>Product weight (without accessories)</b>	100 kg	108 kg
<b>Width</b>	920 mm	1060 mm
<b>Height (floor to top of mattress support platform)</b>	220 mm to 750 mm	220 mm to 750 mm
<b>Length</b>	2030 mm	2030 mm
<b>Safe Working Load (SWL)</b>	200 kg (User weight + mattress weight + accessories weight)	200 kg (User weight + mattress weight + accessories weight)
<b>Max User Weight (MUW)</b>	160 kg (Typical)	160 kg (Typical)
<b>Trendelenburg &amp; Reverse Trendelenburg</b>	0 - 14°	0 - 14°
<b>Backrest Angle</b>	0 - 75° ± 3°	0 - 75° ± 3°
<b>Leg Rest Angle</b>	0 - 45° ± 3°	0 - 45° ± 3°
<b>Castors</b>	8 Swivel Locking Castors (75 mm / 3")	8 Swivel Locking Castors (75 mm / 3")
<b>Input Voltage</b>	100 - 240V	100 - 240V
<b>Output Voltage</b>	DC 29V	DC 29V
<b>Current</b>	3.5A	3.5A

**Specification**
**Power Consumption**
**Long Single**

50W

**King Single**

50W

**Available Attachments**

The Kcare Alivio Jannali Homecare Bed combines innovative functionality with user-friendly design, offering a reliable and stylish solution for individuals requiring supportive care at home.