




KCARE BED STICK WITH SAFETY RETURN

\$82.00



SKU: 9080

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Bed stick with 180 degree safety return . Designed to assist people getting in and out of bed . Gives the advantage of multiple hand positions . White powder coated mild steel . Slides under the mattress for easy installation . Can be used on any size bed and on the left or right

FEATURES

- Bed stick with 180 degree safety return.
- Designed to assist people getting in and out of bed.
- Gives the advantage of multiple hand positions.
- White powder coated mild steel.
- Slides under the mattress for easy installation.
- Can be used on any size bed and on the left or right.

SPECIFICATIONS

Overall Width	500 mm
Overall Height	650 mm
Overall Depth	710 mm