




KCARE BLACK FOREARM CRUTCH WITH ADJUSTABLE HANDLE



SKU: KA343F

PRODUCT DESCRIPTION

Purpose & Benefits of the K Care Forearm Crutch All Size - Black

The K Care Forearm Crutch All Size - Black is designed to provide reliable support and comfort for individuals requiring assistance with mobility. Its multi-adjustable frame accommodates various user sizes and preferences, ensuring optimal stability and ease of use. Ideal for daily use, this crutch enhances mobility and independence.

Features

- **Multi-Adjustable Frame Design:** Allows for independent adjustment of handle and cuff support heights, catering to individual user needs.
- **Lightweight Aluminium Frame:** Constructed from durable yet lightweight aluminium, making it easy to handle and manoeuvre.
- **Soft Padded Hand Grip and Arm Pad:** Provides enhanced comfort during use, reducing strain on hands and arms.
- **Non-Slip Foot Tip:** Ensures secure footing on various surfaces, enhancing user safety.
- **Retractable Design:** Facilitates easy storage and transport, making it convenient for users on the go.

Specifications

Handle Width	136 mm
Min Handle Height	620 mm
Max Handle Height	850 mm

Min Overall Height	890 mm
Max Overall Height	1350 mm
Product Weight	0.7 kg
Max User Weight	136 kg

The K Care Forearm Crutch All Size – Black offers a versatile and comfortable solution for individuals seeking dependable mobility support.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au