



KCARE BLACK FOREARM CRUTCH WITH ADJUSTABLE HANDLE

\$93.00



SKU: KA343F

PRODUCT DESCRIPTION

Purpose & Benefits of the K Care Forearm Crutch All Size - Black

The K Care Forearm Crutch All Size - Black is designed to provide reliable support and comfort for individuals requiring assistance with mobility. Its multi-adjustable frame accommodates various user sizes and preferences, ensuring optimal stability and ease of use. Ideal for daily use, this crutch enhances mobility and independence.

Features

- **Multi-Adjustable Frame Design:** Allows for independent adjustment of handle and cuff support heights, catering to individual user needs.
- **Lightweight Aluminium Frame:** Constructed from durable yet lightweight aluminium, making it easy to handle and manoeuvre.
- **Soft Padded Hand Grip and Arm Pad:** Provides enhanced comfort during use, reducing strain on hands and arms.
- **Non-Slip Foot Tip:** Ensures secure footing on various surfaces, enhancing user safety.
- **Retractable Design:** Facilitates easy storage and transport, making it convenient for users on the go.

Specifications

Handle Width	136 mm
Min Handle Height	620 mm
Max Handle Height	850 mm

Min Overall Height	890 mm
Max Overall Height	1350 mm
Product Weight	0.7 kg
Max User Weight	136 kg

The K Care Forearm Crutch All Size – Black offers a versatile and comfortable solution for individuals seeking dependable mobility support.