



KCARE CLOSED FRONT PU COMMODE SEAT - 445MM

\$178.00



SKU: KA11CPU

PRODUCT DESCRIPTION

FEATURES

- Polyurethane injection moulded seat.
- Includes pan carrier runners.
- Cost effective, comfortable economy seat.
- Clip-on design allows quick removal for cleaning and servicing.
- Suits 445mm wide mobile shower commodes only.
- Open and closed front seats available.

SPECIFICATIONS

| | |
|----------------|--------|
| Overall Width | 440 mm |
| Overall Depth | 430 mm |
| Overall Height | 40 mm |
| Aperture Width | 200 mm |
| Product Weight | 3 kg |
| Brand | K Care |

View the [KCare Tilt In Space Shower Commode](#) or our entire range of [shower commodes](#).