



KCARE DAR-AVOID PRESSURE CARE COMMODE SEAT

\$1,125.00 - \$1,331.55 Price range: \$1,125.00 through \$1,331.55



SKU: KA11CPM-O-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The 'DAR-AVOID' pressure care commode seat is an innovative solution for those who require a higher level of pressure relief. This seat is designed to offload pressure from the ischial tuberosity (IT) and coccyx areas, two of the most commonly affected areas for pressure ulcers. The seat is available in either Dartex or vinyl coverings and features gentle contouring to increase surface area and provide even pressure distribution. Additionally, the 'DAR-AVOID' seat is available with a front opening only, making it easy and convenient for personal hygiene needs. With its superior pressure relief and innovative design, the 'DAR-AVOID' pressure care commode seat is an excellent choice for those who require extra support and comfort while using a commode.

FEATURES

- One of the highest levels of pressure relieving commode seats on the market.
- Designed to offload pressure from ischial tuberosity (IT) and coccyx areas.
- Available covered in dartex or vinyl.
- Gentle contouring to increase surface area and offload pressure.
- Pan runners included.

SPECIFICATIONS

| | |
|----------------|-----------------|
| Seat Width | 445 mm - 550 mm |
| Seat Depth | 460 mm |
| Aperture Width | 200 mm |