



KCARE DELUXE SEAT WALKER 6" - BLUE

\$199.00



SKU: KA365BD

PRODUCT DESCRIPTION

Purpose & Benefits

The **Deluxe Seat Walker 6" - Blue** offers lightweight mobility support with a durable aluminium frame and comfortable seating. Designed for easy indoor movement, it features a padded seat and backrest, smooth soft-touch brakes, and a foldable design for compact transport and storage.

Features

- Strong yet lightweight aluminium frame
- Soft-touch braking system for smooth control
- Easy-fold padded backrest for user comfort and convenience
- Soft padded seat provides comfortable resting support
- Includes under-seat storage bag
- Compact folding design for easy transport and storage
- Intended for indoor use only

Specifications

Specification	Detail
Overall Length	650 mm
Overall Width	620 mm
Min Overall Height	840 mm
Max Overall Height	940 mm
Seat Width	445 mm
Floor to Seat Height	550 mm
Castor Diameter	150 mm (6")
Product Weight	6.2 kg
Max User Weight	160 kg
Brand	K Care