


**KCARE DELUXE SEAT WALKER 6" - BLUE****SKU:** KA365BD**PRODUCT DESCRIPTION****Purpose & Benefits**

The **Deluxe Seat Walker 6" - Blue** offers lightweight mobility support with a durable aluminium frame and comfortable seating. Designed for easy indoor movement, it features a padded seat and backrest, smooth soft-touch brakes, and a foldable design for compact transport and storage.

**Features**

- Strong yet lightweight aluminium frame
- Soft-touch braking system for smooth control
- Easy-fold padded backrest for user comfort and convenience
- Soft padded seat provides comfortable resting support
- Includes under-seat storage bag
- Compact folding design for easy transport and storage
- Intended for indoor use only

## Specifications

Specification	Detail
Overall Length	650 mm
Overall Width	620 mm
Min Overall Height	840 mm
Max Overall Height	940 mm
Seat Width	445 mm
Floor to Seat Height	550 mm
Castor Diameter	150 mm (6")
Product Weight	6.2 kg
Max User Weight	160 kg
Brand	K Care

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
 Brighton SA 5048  
 08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
 Ave, Elizabeth Vale SA 5112  
 08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
 Reynella SA 5161  
 08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
 03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
 03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)