



KCARE ECONOMY BEDSIDE COMMODE

\$238.00



SKU: KA500PE03-H

PRODUCT DESCRIPTION

Purpose & Benefits of the KCare Economy Bedside Commode

The KCare Economy Bedside Commode is designed to provide a convenient and comfortable toileting solution for individuals with limited mobility. It allows for safe and easy access to toileting facilities, especially during the night or in environments where traditional bathrooms are not easily accessible.

Features

- **Corrosion-Resistant Frame** – Constructed from phosphate-treated mild steel with an epoxy powder coating, ensuring durability and longevity.
- **Padded Support** – Equipped with a padded vinyl backrest and removable padded seat, offering enhanced comfort during use.
- **Adjustable Height** – Features height-adjustable legs to accommodate various user needs and preferences.
- **Supportive Armrests** – Provides additional support during sitting and transferring, enhancing user safety.
- **Removable Components** – Includes a lift-up toilet seat and a fixed bowl carrier for easy cleaning and maintenance.
- **Autoclave Bowl** – Supplied with a bowl that can be autoclaved up to 120 degrees Celsius, ensuring hygiene and ease of cleaning.

Specifications

- **Overall Width:** 560 mm
- **Overall Depth:** 590 mm
- **Seat Width:** 440 mm

- **Seat Depth:** 440 mm
- **Minimum Seat Height:** 490 mm
- **Maximum Seat Height:** 615 mm
- **Backrest Height:** 450 mm
- **Product Weight:** 6.5 kg
- **Maximum User Weight:** 110 kg

The KCare Economy Bedside Commode is an ideal solution for individuals seeking a reliable and comfortable bedside toileting option, combining practicality with user-friendly features.