



KCARE ECONOMY BEDSIDE COMMODORE



SKU: KA500PE03-H

PRODUCT DESCRIPTION

Purpose & Benefits of the KCare Economy Bedside Commode

The KCare Economy Bedside Commode is designed to provide a convenient and comfortable toileting solution for individuals with limited mobility. It allows for safe and easy access to toileting facilities, especially during the night or in environments where traditional bathrooms are not easily accessible.

Features

- **Corrosion-Resistant Frame** – Constructed from phosphate-treated mild steel with an epoxy powder coating, ensuring durability and longevity.
- **Padded Support** – Equipped with a padded vinyl backrest and removable padded seat, offering enhanced comfort during use.
- **Adjustable Height** – Features height-adjustable legs to accommodate various user needs and preferences.
- **Supportive Armrests** – Provides additional support during sitting and transferring, enhancing user safety.
- **Removable Components** – Includes a lift-up toilet seat and a fixed bowl carrier for easy cleaning and maintenance.
- **Autoclave Bowl** – Supplied with a bowl that can be autoclaved up to 120 degrees Celsius, ensuring hygiene and ease of cleaning.

Specifications

- **Overall Width:** 560 mm
- **Overall Depth:** 590 mm
- **Seat Width:** 440 mm

- **Seat Depth:** 440 mm
- **Minimum Seat Height:** 490 mm
- **Maximum Seat Height:** 615 mm
- **Backrest Height:** 450 mm
- **Product Weight:** 6.5 kg
- **Maximum User Weight:** 110 kg

The KCare Economy Bedside Commode is an ideal solution for individuals seeking a reliable and comfortable bedside toileting option, combining practicality with user-friendly features.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au