



KCARE ECONOMY PADDED STEEL SHOWER STOOL HEIGHT ADJUSTABLE

~~\$100.00~~ Original price was: \$100.00. ~~\$90.00~~ Current price is: \$90.00.



SKU: 12325S

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Economy Shower Stool with Padded Seat** provides a safe and comfortable solution for users who experience difficulty standing while showering. Its padded PU seat offers extra comfort, while the durable steel frame ensures strength and stability. This stool promotes independence and reduces the risk of slips and falls in the bathroom.

FEATURES

- Assists individuals who have difficulty showering independently
- Soft PU padded seat provides enhanced comfort and support
- Height-adjustable legs to suit a wide range of users
- Powder-coated 22mm steel tubing for durability and corrosion resistance
- Compact design ideal for smaller bathrooms
- Lightweight and easy to reposition when needed

SPECIFICATIONS

- **Min Overall Height:** 680 mm
- **Max Overall Height:** 810 mm
- **Seat Width:** 410 mm
- **Seat Depth:** 300 mm
- **Seat Height Range:** 460 – 600 mm

- **Max User Weight:** 110 kg
- **Frame Material:** Powder-coated 22mm steel tubing
- **Seat Material:** PU padded