


KCARE ECONOMY SHOWER CHAIR



SKU: 12028S



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **KCare Economy Shower Chair** provides a safe and stable seating option for users who require additional support while showering. Its strong steel frame and supportive armrests assist with positioning and transfers, helping to reduce the risk of slips and falls in the bathroom.

FEATURES

- Provides secure seating while showering
- Armrests offer extra stability and assistance with sitting and standing
- Constructed from durable 22mm powder-coated white steel tubing
- Solid plastic seat for strength and easy cleaning
- Height adjustable legs for personalised comfort
- Designed and tested to Australian standards for safety and reliability

SPECIFICATIONS

Specification

Detail

Min Overall Height	610 mm
Max Overall Height	730 mm
Seat Width	460 mm

Specification

Detail

Seat Depth	480 mm
Min Seat Height	430 mm
Max Seat Height	550 mm
Maximum User Weight	110 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au